



Select One Service Level	Academic Service <input type="checkbox"/>	Calibration <input type="checkbox"/>	Gold* Renewal <input type="checkbox"/>
Single Channel Pipettes	\$26	\$35	\$85
Electronic Single Pipettes	\$30	\$40	\$125
Repeater Pipettes (manual/electronic)	\$30	\$40	\$85/\$155
8-Channel Pipettes (manual/electronic)**	\$45	\$60	\$155/\$215
12-Channel Pipettes (manual/electronic)**	\$55	\$70	\$155/\$255
Motorized Controller	\$35	\$40	N/A

Service Includes			
Leakage Test & Basic Cleaning	included	included	included
Functional & Calibration Check	included	included	included
Preventative Maintenance (PM) ***	included	included	included
All Replacement Parts, as needed	extra	extra	included
CAL Weigh Readings	2x2	4x2	4x2
Service Report	Academic	Cal Report	GMP/GLP
Comprehensive CAL Label	included	included	included
Free Shipping for this minimum quantity	ten	five	five

Optional Service Features			
Printed Copy of GMP/GLP Report****	n/a	<input type="checkbox"/> \$3	<input type="checkbox"/> JbWl XYX
DNA/RNA Decontamination "per pipette"	<input type="checkbox"/> \$10	<input type="checkbox"/> \$10	<input type="checkbox"/> %\$

* Pipette.com reserves the right to exclude some pipettes from the Gold Renewal program
 ** Calibration on 8 & 12 channel pipettes is based on an average weightings of all channels.
 *** PM includes new seal and o-ring as specified by manufacturer, piston polishing, thorough inside and outside cleaning, full refurbishment, no labor charge for parts installation.
 **** GMP/GLP reports are available on line at no charge.
 All pipette brands serviced and repaired in our environmentally controlled ISO 17025 Accredited Lab using NIST ultra-micro (up to 10⁻⁷g) traceable balances.
 Calibration Warranty: 90 days labor / 1 year installed parts.

For SHIPPING instructions go to www.PIPETTE.COM/SHIPPING.

Shipping & Handling Options

- Request 48 hours service Free
- 24 Hour Service (per pipette) \$5
- Priority Overnight Shipping (each way) \$28

List pipettes below or attach your pipette list to this form.

Model Number	Serial Number	User/ID/Notes

QUANTITY OF PIPETTES SHIPPED: NEXT CALIBRATION DUE IN Months

Pipettes that require rare parts may be delayed in returned shipment.

FULL NAME

PHONE

E-MAIL

FAX

SHIP TO ADDRESS

COMPANY NAME

DEPARTMENT

ADDRESS

CITY, STATE & ZIP CODE

BILL TO ADDRESS

COMPANY NAME

DEPARTMENT

ADDRESS

CITY, STATE & ZIP CODE

PAYMENT INFORMATION

- MASTER CARD VISA AMERICAN EXPRESS

PO NUMBER OR CREDIT CARD NUMBER

EXP DATE

PIPETTE.COM ACCOUNT (optional)

SERVICE CODE (optional)

BILL CUSTOMER SHIPPING ACCOUNT 'list info below'

Special Requests

PIPETTES WERE DECONTAMINATED FROM RADIOACTIVE, CHEMICAL AND BIOHAZARDOUS MATERIALS WITH AN APPROPRIATE SOLUTION

CUSTOMER SIGNATURE